VOCO

Alfacomp

(corresponde a EN 24049/ISO 4049)

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(in accordance with EN 24049/ISO 4049)

Instructions for Use

Alfacomp is a BIS-GMA-based, chemical curing composite with 81% w/w filler content (ca. 0.05 and 3 - 5 μ). **Alfacomp** comes in a universal shade (A3) which is suitable for most patients.

Indications:**Alfacomp** is recommended for:

- class III - V fillings
- class I and II fillings which are not exposed to occlusal forces
- extended fissure sealing
- core build-up
- reconstruction of traumatically affected teeth
- improving aesthetics (e.g. diastema, hypoplasia, conical teeth, erosions)

Application:

Clean the teeth to be treated with a fluoride-free cleaning paste. Using a minimally invasive technique/adhesive filling technique is recommended to preserve healthy tooth substance.

Cleaning and drying the cavity has to be carried out with special preparations (e.g. hydrogen peroxide solution). Observe instructions for use when applying the bonding materials.

Apply pulp protection when working in close proximity to the pulp. Lining deep cavities with a glass ionomer cement is recommended due to its release of fluoride. Bevel enamel margins and etch for 20 - 30 s (**Vococid**). Then thoroughly rinse and dry. Do not etch dentine for more than 15 s when using the "total-etch-technique". Use a suitable dental adhesive according to the instructions for use.

Mix equal quantities of **Alfacomp** base and catalyst paste on the mixing pad with a disposable spatula for approx. 15 - 30 s until homogeneous.

Do not use the same end of the spatula for base and catalyst paste. If any traces of one paste get into the other, the composite will become useless. Dry the cavity before application.

The working time of **Alfacomp** is approx. 2 min and the setting time is 3 min. Do not move strips or matrix during this time.

Remove any large areas of excess material with a diamond bur. Using discs of graduated grains or a finishing diamond (lingual, palatinal area) is recommended for finishing. Diamond metal strips are suited for the approximal area. Then polish with a suitable polishing paste.

A fluoride varnish should be applied to all treated teeth to support remineralization.

Warnings and precautions:

- In case of contamination, etch the already etched surfaces again for 5 s, rinse and dry.
- Do not use light-curing one-bottle bonds with **Alfacomp**.
- Alfacomp** contains BIS-GMA, BHT, benzoyl peroxide; some patients (e.g. allergies) may be sensitive to these ingredients.
- Alfacomp** should not be used if the cavity cannot be dried, the prescribed technique cannot be used, or if there are known allergies to its ingredients.
- Preparations containing eugenol interfere with the setting of composite restoratives. Therefore, the use of zinc-oxide eugenol cements or other preparations containing eugenol should be avoided.
- Cationic mouthwashes, plaque detectors and chlorhexidine may discolour the restoration and should be avoided.

Storage:

Store at temperatures between 4°C - 23°C. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

Last revised: 10/2009	
VOCO GmbH P.O. Box 767 27457 Cuxhaven Germany	Phone +49 (4721) 719-0 Fax +49 (4721) 719-140 e-mail: marketing@voco.com www.voco.com

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VC 60 001070 E1 109 99 © by VOCO