VOCO

Aqua Ionofil Plus

in accordance with **GE** **DIN EN 29917** **CE** ⁰⁴⁸² **Instructions for use**

Aqua Ionofil Plus is a radiopaque, water-mixable glass ionomer filling cement with chemical adhesion to enamel and dentine. **Aqua Ionofil Plus** is available in shades A1, A3 and B3 and has enamel-like translucency.

Indications:

- for restoration of deciduous teeth (especially class I)
- repair of V-shaped defects and enamel as well as root erosions in the dental neck area including class V fillings
- cavity fillings with carious lesions
- fillings of class III and small fillings of class I (incl. extended fissure sealing) as well as core build-up under crowns and as a cavity lining

Method of Use:

The adhesion of **Aqua Ionofil Plus** to enamel and dentine prevents the formation of marginal leakage and creates retentions. Undercuts provide additional hold to the filling. Layer thickness of the filling should not be less than 1 mm.

Clean the cavity (e.g. with a 3 % hydrogen peroxide), rinse carefully with water and dry. In case of remaining dentine thickness of more than 1 - 1.5 mm no additional lining is necessary. In case of inferior dentine thickness it is recommended to cover the areas in proximity of the pulp with a calcium hydroxide.

Attention! Only mix Aqua Ionofil Plus with water.

Shake powder well before use.

Allow **Aqua Ionofil Plus** to reach room temperature (max. 23°C) before application. Use a cooled glass plate if necessary. Dispense **Aqua Ionofil Plus** powder on a glass plate or suitable mixing pad and mix the powder portion by portion with the water. Use 1 drop of water for 1 level measuring scoop of powder (mixing ratio powder : water = 5.4 - 6.3 g : 1).

When dispensing the water, hold the dropper bottle in a vertical position, making sure that the metal cannula protrudes approx. 2 - 3 mm from the dropper. Do not use drops for mixing which escape when bending the bottle. Use a strong plastic spatula or another abrasion-resistant spatula for mixing.

| | | |
|--|--|--|
| Mixing time | Working time | Setting time |
|  |  |  |
| 30 - 40 s | 2.5 - 3.5 min | 5 - 6 min |

The mixing time is approx. 30 - 40 s, but should not exceed 1 min.

Working time from end of mixing procedure is approx. 2 - 3 min (at room temperature). After the filling has been inserted into the cavity, the setting time of **Aqua Ionofil Plus** is approx. 2 - 3 min.

Ensure that the cavity remains dry during the filling procedure.

Class V cavities should be at least 1 mm in depth. Prevent feather-shaped margins. In case of doubt prepare in a conventional form.

During the setting phase, especially immediately after application, glass ionomer cements are sensitive to humidity (breath, saliva) and drying-out. Therefore, cover immediately with a strip or matrix or apply a protective varnish (see instructions for use **Final Varnish LC**) if possible.

The rough excess can be removed with a sharp instrument after **Aqua Ionofil Plus** has set (4 - 6 min after application depending on mixing temperature and ratio). Finish with a fine-grained diamond or a flexible disk. Afterwards, apply a coating of varnish. For best results, perform final polishing 24 h after placing the filling.

Precautions:

In order to avoid irritations in the proximity of the pulp (less than 1 - 1.5 mm) use pulp protection. Tightly close powder after use and protect from moisture.

Avoid contact with the mucosa and eyes.

Aqua Ionofil Plus contains fluorosilicate and polyacrylic acid. Do not use in case of sensitivities to these ingredients, after expiry date or when the powder forms lumps. Do not use insulators on dentine in order to provide maximum adhesion with the exception of pulp protection.

Storage:

Store at temperatures between 4°C - 23°C.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VOCO GmbH **Phone +49 (4721) 719-0**