VOCO**Arabesk**

in accordance with EN 24049/ISO 4049

Instructions for use

Arabesk is a universal, light-curing, radiopaque micro hybrid dental composite filling material for the anterior and posterior areas, as well as for the inlay technique. It contains 60 % by volume (= 76.5 % by weight) inorganic fillers, microfillers (approx. 0.05 µm) and small particle fillers (approx. 0.5 - 2 µm). **Arabesk** can be polished to a high gloss and is distinguished by its high level of stability and colour stability. **Arabesk** is polymerised by halogen light (blue light) or LED-light.

Arabesk is available in practical syringes for direct application.

It is available in 7 shades:

reddish/brownish shades	A1, A2, A3, A3.5
yellowish shades	B2, B3
incisal	I

Fields of application:

- fillings of class I, II, V in the posterior area which are not exposed to occlusal forces
- reconstruction of traumatically damaged anteriors
- facing of discoloured anteriors
- shape and shade corrections to improve aesthetic appearance
- locking of loosened anteriors
- extended fissure sealing for molars and premolars
- fillings of classes III, IV, V in the anterior area
- repair of veneers
- core build-ups for crowns
- especially suitable for composite inlays

Methods of use:

Preparation/shade selection: Clean the teeth to be treated with fluoride-free cleaning paste. Mark occlusal contact points (posterior area), a minimal separation makes the shaping of the approximal contacts and the placement of the matrix easier. For shade matching compare the **Arabesk** shade guide made of original material (corresponds to the shade of the filling after placement and sufficient polymerisation) with the cleaned tooth before anaesthesia (moisten the shade guide and the tooth and view in daylight). Use of a rubber dam is recommended.

Cavity preparation:

In principle, minimal preparation (adhesive filling technique) should be carried out, preserving the healthy tooth substance. Special preparation methods have to be observed for extended fissure sealing and for primary teeth!

Cleaning/Drying:

The cleaning and drying of the cavity has to be carried out with special preparations or hydrogene peroxide solution. When applying bonding materials observe their specific characteristics (see corresponding instructions for use).

Lining:

A calcium hydroxide lining should be applied in proximity of the pulp. In case of very deep cavities, lining is recommended. Due to their fluoride release, glass polyalkenoate (glass ionomer) cements are most suitable for this purpose.

Etching/Bonding:

Use **Arabesk** with the "total-etch-technique" and a bonding agent (e.g. **Solobond M**). Please follow the instructions for use of **Vococid** and **Solobond M**.

Application of Arabesk:

Allow the material to reach room temperature before application. Select the matching shade of **Arabesk** (see shade matching) and apply with a suitable instrument. It is advantageous to use translucent matrices (smooth surfaces).

Further treatment:

For fillings of more than 2 mm in thickness, apply and polymerise layer by layer (cure every layer for approx. 40 - 60 s). To light-cure this material conventional polymerisation devices are appropriate. The light output should be a minimum of 500 mW/cm² on halogen polymerisation devices and 300 mW/cm² on LED devices. Transparent strips, crowns or light wedges should be used. Care should be taken to ensure good marginal sealing. Hold the light source of the halogen light as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth will be smaller. It is then expected that the polymerisation will be limited to the area of the cone of light only. Insufficient curing can lead to discolouration and pulp irritation.

The filling can be finished and polished immediately after removal of the transparent strips/crowns (fine or extra-fine grain diamond, polishing discs). During finishing and polishing the filling must be cooled. Finally, the margin of the filling, or preferably the entire tooth, should be fluoridated.

Arabesk can be used for direct and indirect inlays according to the known methods. The physical stability values can be improved by the usual external (additional) curing.

Information:

- No known side-effects. However, sensitisation of hyper-sensitive persons cannot be ruled out.
- If the cavity cannot be dried or the prescribed technique cannot be used, or in case of known allergies to certain ingredients (Bis-GMA, UDMA, TEGDMA, BHT, HEMA), **Arabesk** should not be used.
- Materials containing eugenol, thymol or oil of cloves inhibit the setting of filling composites. The use of zinc oxide eugenol cements or other materials containing eugenol in combination with filling composites should therefore be avoided.
- When in contact with cationic mouthwashes, plaque solution or chlorhexidine, discolouration may occur.

Storage:

Store at temperatures of 4°C to 23°C. Syringes and bonding material bottles should be closed immediately after material has been extracted to prevent exposure to light. Do not use after the expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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