VOCO

Argion Molar AC

GB **CE** 0482 **in accordance with EN 29917/ISO 9917**

Instructions for use

Product description:

Argion Molar AC is a silver-containing, radiopaque glass polyalkenoate (ionomer) cement in Application Capsules (AC) for direct application.

Argion Molar AC chemically adheres to enamel and dentine and continuously releases fluoride ions.

0.15 ml cement can be extracted from each capsule.

Indications:

- restoration of deciduous teeth (especially class I)
- class II filling (except occlusion-bearing) by means of tunnel preparation
- build-up fillings, linings, core build-up
- long-term provisional treatment of class I and II cavities


Type of application:

Preparation:

The adhesion of **Argion Molar AC** to enamel and dentine prevents the formation of marginal gaps and creates retentions.

Undercuts give additional hold to the filling. The layer thickness of fillings should not be less than 1 mm. Avoid feather-shaped margins. If in doubt, choose conventional preparation.

Remove preparation residues with water spray, then dry with air so that the surface remains slightly glossy (do not overdry). If the remaining dentine is thicker than 1.5 mm, additional lining is not necessary. Cover areas in close proximity to the pulp with calcium hydroxide to avoid irritations which may occur in cases of insufficient remaining dentine thickness. The capsule must be activated, mixed and opened before application.

Activation:

For activation, the liquid from the sachet underneath the clamp has to be emptied as completely as possible into the inner chamber of the capsule.

This is done with a capsule activator. Place the capsule in the activator so that the application cannula protrudes from it and the clamp faces the moveable part of the activator. **Hold the lever for 2 s in completely compressed position.**

Mixing:

Fix the activated capsule in a high-frequency capsule mixer (4,000 - 5,000 oscillations per min) and mix for approximately 10 s. Immediately remove the capsule from the mixer.

The use of capsule mixers of other manufacturers may result in shorter or longer mixing times due to different mixing frequencies.

Opening:

Open the twist cap with application cannula with a half turn to the right until you notice a stop (see stamp on protective cap).

Application:

Insert the opened capsule into the holding bracket of an applicator. By turning the capsule in the holding bracket the cannula can be aligned in the desired direction.

Apply the material within 30 s of mixing it directly into the cavity by pressing the lever of the applicator.

Mixing time	working time	setting time
0 - 10 - 15 s	1.5 - 2.5 min	3.5 - 4.5 min

After mixing, **Argion Molar AC** has a working time of approximately 1.5 min (room temperature). After application in the cavity **Argion Molar AC** sets within 2 - 3 min. During the placing of a filling the work field must be kept dry. The use of a rubber dam is recommended. Linings and build-up fillings can be finished carefully when set (after approximately 4 - 6 min in the mouth). Glass ionomer cements are susceptible to humidity (breath, saliva) or, respectively, dehydration during the setting time, especially directly after application. If possible, cover immediately with a strip or matrix or apply a protective varnish (see instructions for use **Final Varnish LC**). Finishing and polishing (flexible polishing discs) of fillings should be conducted no earlier than 7 min after placement. Then re-apply protective varnish. The best results of fine finishing will be achieved approximately 24 h after placement of the filling.

Additional advice/precautions:

Argion Molar AC should be mixed only immediately before application in order to have the full working time available for inserting and shaping and to avoid cracks from working the material too late.

Please note that the AC capsule has to be mixed right after its activation.

Argion Molar AC contains polyacrylic acid, fluoroalicates, silver and parabenes. Do not use in cases of known sensitivities or allergies.

Storage:

Store **Argion Molar AC** mixing capsules in a dry area at temperatures between 15 °C - 23 °C. Do not remove foil until ready to use. Do not use after the expiry date printed on this package. Tightly close **Final Varnish LC** container immediately after use. Store at room temperature and protected from light.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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