VOCO

Calcimol LC

(GB)

Instructions for use

Calcimol LC is a light-curing radiopaque one-component material containing calcium hydroxide.

Indications:

- Indirect pulp capping
- As a liner under all filling materials
- As a protection against acid when using the etch technique

Application:

The dental surfaces have to be dried thoroughly and contamination avoided.

The preparation complies with the rules of the adhesive restorative technique. Use a diluted calcium hydroxide preparation instead of **Calcimol LC** in direct proximity to the pulp.

To prevent marginal leakage and infiltration of bacteria in the cervical/approximal margins, permanent isolation is recommended.

Calcimol LC can be directly applied from the syringe.

In order to ensure the function of the non-dripping syringe do not pull back the plunger during or after use.

From tubes, dispense the required quantity of **Calcimol LC** onto a mixing pad and apply with a suitable instrument. In case of linings of more than 1 mm apply and polymerise layer by layer.

Cure **Calcimol LC** with blue light. Keep the light source as close as possible to the surface (approx. 1 - 2 mm) and cure for at least 20 s. In case of more distance extend curing time. Carefully remove excess material with a rotating instrument and apply the filling according to the manufacturer's instructions.

Warnings/precautions:

Calcimol LC adheres to materials, e.g. crown and bridge materials, temporary sealants on methacrylate basis, certain impression materials (e.g. polyether). This could lead to removal of the lining (core build-up). In such cases it is recommended to apply an insulating agent (e.g. vaseline, glycerine gel or similar). Remove insulating agent prior to definite application.

The consistency of **Calcimol LC** does not allow to remove excess material with a probe, this could eventually remove the whole lining (use a rotating instrument).

After use of **Calcimol LC** avoid direct light (operation lamp, sun-light) and apply the material immediately.

Calcimol LC contains calcium hydroxide, TEGDMA, BHT, UDMA and amine. Do not use in patients with known allergies to these ingredients.

Eugenol and other contaminates (e.g. thymol) prevent the polymerization and therefore should be avoided. Insufficient curing times and/or light intensity lead to incomplete curing. Therefore the halogen lamp/light guide should be checked regularly. Tightly close after use and store in an area protected from light and moisture at temperatures between 4°C - 23°C.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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