VOCO Final Varnish LC

GB Instructions for use

Product description:
Final Varnish LC seals and smoothes the surface of glass ionomer cements. The light-curing covering varnish **Final Varnish LC** produces a smooth surface and protects the cement surface from damage due to drying or ingress of moisture or saliva until the cement is cured completely.

Indications:

- protection of glass ionomer cement surfaces from the effects of humidity and dryness immediately after placing the filling
- sealing and smoothing of glass ionomer cement surfaces after polishing
- sealing of the cervical margins of cementations immediately after removal of excess material
- sealing of linings/build-up fillings before taking impressions

Directions:
Dispense 1 - 2 drops onto a pad or application instrument and close bottle immediately.

Apply **Final Varnish LC** with a brush or another appropriate instrument immediately after placing the filling or finishing the glass ionomer cement. Remove any air bubbles with a probe. Cover all glass ionomer surfaces with the varnish (approximal surfaces or surfaces after removal of matrices should also be covered). Cure **Final Varnish LC** with halogen light/LED-light (10 s per surface). To light-cure this material conventional polymerisation devices are suited. The light output should be a minimum of 500 mW/cm² for halogen polymerisation and 300 mW/cm² for LED devices.

If the varnished area requires subsequent treatment within 24 h after the initial application, reseal it afterwards.

Precautionary measures/advice Final Varnish LC:

Final Varnish LC contains Bis-GMA, diurethanedimethacrylate, BHT.

Do not use on patients with known hypersensitivities/allergies to any ingredients.

To avoid its polymerisation (increased viscosity), expose **Final Varnish LC** to the chairside light for a short time only. Check sulcus after application, remove any excess material.

Avoid contact with mucosa.

Apply **Final Varnish LC** again, if the seal is damaged prior to final polymerisation.

Unpolymerised or incompletely dried varnish may impede setting of addition-curing silicones.

Storage:

Tightly close **Final Varnish LC** bottle immediately after use. Store at room temperature and protected from light. Do not use after expiry date.

Final Varnish LC is a registered trademark of VOCO. Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 AA1241 E1 1211 99 © by VOCO