**VOCO**

# Futurabond NR

**GB**

Instructions for use

**CE** 0482

**Futurabond NR** is a light curing self-etch-bond reinforced with nanoparticles to build a durable bond between tooth substance and light-curing filling materials without marginal leakage. This self-etching adhesive consists of 2 components. A separate etching step is not necessary: **Futurabond NR** shows the same adhesive properties as total-etch bonds. **Futurabond NR** tolerates residual moisture.

**Futurabond NR** can be used universally with all light-curing composites, compomers and Ormocer<sup>®</sup>s.

**Indications:**

**Futurabond NR** is suited for direct restorations with light-curing filling materials for all cavity classes.

**Directions:**

**1. Drying**

Due to the subsequent application of light-curing filling materials a dry work field is required. Avoid any contamination with blood or saliva.

**2. Preparation**

Prepare the cavity according to the principles of adhesive restoration technique. Rinse cavity thoroughly with water. Remove excess moisture with faint air jet. Do not overdry dentine. The dentinal surface should remain slightly moist. If the tooth has not been prepared (e.g. cervical defect) the surface has to be cleaned thoroughly and the enamel margin roughened with a finishing diamond.

Areas close to the pulp should be protected with suitable calcium hydroxide-based liner.

**3. Mixing**

**Futurabond NR** is light-curing and therefore has to be applied immediately after mixing.

**When using the bottle:**

Mix 1 drop of Liquid **A** and 1 drop of Liquid **B** on a mixing palette with a **Single Tim** for approximately 5 s (produces a self-etching adhesive)

**When using *SingleDose*:**

Detach a *SingleDose* by folding along the perforation and pulling the printed side upwards. To activate, hold the *SingleDose* blister between thumb and forefinger and press on the blister where indicated “press here” in blue. This injects Liquid **B** into Liquid **A** in the dispensing bay in a 1:1 ratio. Pierce the foil of the dispensing bay with the **Single Tim**. Enlarge the hole with a circular movement and create a homogeneous mixture of A and B by repeated stirring.

**4. Application of Futurabond NR**

Apply a moderately thin layer of the adhesive to the enamel/dentine and massage into the tooth substance for 20 s.

**5. Drying**

Dry the adhesive layer with an air jet for at least 5 s.

**6. Polymerization**

Polymerize with blue light (halogen- or LED-light) for 10 s.

**Advice and precautions:**

- Do not fold blister before activation. Blisters with obvious damage such as creases or kinks should not be used.

- Do not damage either the blister foil or tab on either the top or the bottom side with fingernails.

- Do not use any type of instrument for activation. Only the thumb and index finger should be used for activation with the thumb placed exactly on the label's mark and the index finger placed exactly on the opposite side of the package.

- Activation and subsequent puncturing of the foil with a **Single Tim** should not take place within close proximity of the eyes.

- Never simultaneously press the both reservoirs or press the liquid back and forth between the reservoirs.

- Phenolic substances (e. g. eugenol-containing temporary materials) disturb the polymerization and must not be brought in contact with **Futurabond NR**.

- **Futurabond NR** contains organic acids, Bis-GMA, HEMA, TMPTMA, campherquinone, amines (DABE), BHT, fluorides and ethanol.

Do not use in cases of known allergies to any of these ingredients.

Avoid contact with the oral mucosa. In case of contact of **Futurabond NR** with the eyes, rinse immediately with plenty of water and consult a physician.

**Single Tim** fine application brushes allow a precise and thin application.

**Information on storage and application:**

Tightly close **Futurabond NR** after use. Store in upright position in order to guarantee flow-back of the liquid. Store at temperatures between 4°C - 23°C. Do not expose **Futurabond NR** to direct sunlight or chair-light.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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**VOCO GmbH**                   **Phone +49 (4721) 719-0**