VOCO

Glasiosite

entspricht EN 24049/ISO 4049

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in accordance with EN 24049/ISO 4049

Instructions for use

Glasiosite is a universal light-curing radiopaque compomer restorative material for filling cavities in both the anterior and posterior areas. It allows a simplified and rapid filling technique as well as combines the fluoride release of glass ionomers with the physical durability and high aesthetics of composites. **Glasiosite** is used in combination with **Futurabond NR**.

Futurabond NR is a universal self-etching bond which provides uniformly high adhesion to enamel and dentine.

Indications:

- filling deciduous teeth
- class V fillings (cervical caries, root erosions, V-shaped defects)
- faceting of discoloured anteriors
- class III fillings
- temporary class I and II fillings
- extended fissure sealing

Contraindications:

- **Glasiosite** is contraindicated:
- If there known allergies to its ingredients (methacrylates, amines)
- If a dry work field cannot be secured.
- If the prescribed application technique cannot be employed.
- For direct and indirect pulp capping
- For load-bearing permanent fillings in secondary dentition.

Application:

Preparation/ shade selection:

Clean the teeth to be treated with fluoride-free cleaning paste. To select the shade, compare the shade guide with original material (corresponds to the shade of an applied and thoroughly polymerized filling) with the cleaned tooth before anaesthesia (moisten shade guide and the tooth, select the shade in daylight). Secure a dry working field.

Cavity preparation:

In general, preparation should be minimal (adhesive filling technique) to conserve healthy tooth substance.

Observe special preparation forms for extended fissure sealing and on deciduous teeth.

Caries-free cervical defects do not require preparation. Clean the defect with either cleaning paste or pumice using a rubber cup or rotating brush.

Cleaning:

- Remove all residues from the cavity with the water jet.
- Avoid any contamination with blood or saliva after the cleaning.
- Remove excess water with a faint air jet, yet do not overdry the dentine. A moist dentine surface is necessary for wet-bonding.

Pulp protection:

A calcium hydroxide lining should be applied in close proximity to the pulp.

Application of Futurabond NR:

Futurabond NR is light-curing and therefore has to be applied immediately after mixing.

When using the bottle:

Mix 1 drop of Liquid **A** and 1 drop of Liquid **B** on a mixing palette with a **Micro Tim/ Single Tim** for approximately 5 s (produces a self-etching adhesive)

- Application of Futurabond NR

Apply a moderately thin layer of the adhesive to the enamel/dentine and massage into the tooth substance for 20 s.

- Drying

Dry the adhesive layer with the air jet for at least 5 s.

- Polymerization

Polymerize with blue light (halogen- or LED-light) for 10 s.

Application of Glasiosite:

Apply and light-cure fillings more than 2 mm deep in layers. Polymerize each layer for 40 s with a conventional polymerization device. The light output should be a minimum of 500 mW/cm² on halogen polymerization devices and 300 mW/cm² on LED devices.

Use transparent strips and crowns or light wedges. Ensure a tight marginal seal. Hold the tip of the light rod as close as possible to the surface of the filling. If the distance exceeds 5 mm the curing depth may be reduced to polymerization in the range of the light cone only. Incomplete polymerization may lead to discoloration and pulpitis-like complaints.

Finishing and polishing of the filling can be conducted with cooling provided (e. g. fine or extra fine finishing diamonds, polishing discs) immediately after removing the shaping aids. As the last step, the margin of the filling or the entire tooth should be fluoridated.

Indications:

- **Glasiosite** contains BIS-GMA, diurethandimethacrylate, TEGDMA, BHT. Do not use in patients with known allergies to these ingredients.
- Preparations containing eugenol or oil of cloves disturb the polymerization of resin-based filling materials and should not be used in combination with resin-based filling materials.
- Contact with cationic mouth rinses as well as plaque indicators and chlorhexidine may lead to discoloration.

Storage:

Store at temperatures between 4°C - 23°C. Close syringes and bond bottles immediately after dispensing to avoid premature curing and/or evaporation of the solvent.

Presentations:

Glasiosite Caps

Set 40 x 0.25 g Caps, shades
10 x A3, 10 x A4, 5 x A2, 5 x A3.5,
5 x B3, 5 x C2, 2 x 4 ml **Futurabond NR**
light-curing self-etch-bond
and accessories Art. No. 1600

Refills 25 x 0.25 g Caps, shades
A2, A3, A3.5, A4, B3, C2 Art. No. 1601 - 1606

Shade guide Art. No. 2101

Caps Dispenser Art. No. 9301

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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© by VOCO