VOCO

Solobond Plus

(GB)

Instructions for use

(€)⁰⁴⁸²

Solobond Plus is a light-curing two-stage universal bonding agent for permanent bonding. The **Solobond Plus** system consists of a **primer** to condition the dentine and an **adhesive** as bonding agent between conditioned dentine, etched enamel and composite. **Solobond Plus** yields tight marginal seals for the prevention of marginal caries.

Solobond Plus tolerates residual humidity and permits "Wet Bonding". The special composition of **Solobond Plus** enables the self-conditioning priming (*without* etching the dentine, see method of use A), as well as the use of the total-etch-technique (*with* etching the dentine, see method of use B). Thus shallow cavities can be restored by using the practical total-etch-technique. For deep cavities in proximity of the pulp, the self-conditioning primer is recommended. **Solobond Plus** can be used with both light-curing and self-curing composite materials.

Indications:

Solobond Plus is suitable for direct composite restorations of all classes as well as for adhesive luting of crowns, bridges, inlays, onlays, veneers etc. with composite cements.

Method of use:

A) Direct composite restoration with self-conditioning priming

1. Drying

The application of **Solobond Plus** requires dry tooth surfaces due to the later application of the composite. Avoid contamination with blood or saliva. The use of rubber dam is recommended.

2. Preparation

Prepare the cavity according to the principles of the adhesive filling technique. Do not remove the smear layer! Protect areas in proximity of the pulp with a suitable calcium hydroxide based liner.

3. Etching of enamel margins

Apply **Vococid** etching gel to the areas to be etched and let act for ca. 20 - 30 s. Then remove gel by suction and carefully rinse with water. Remove residual humidity (air syringe).

4. Application of the primer

Dose primer onto a mixing pallet and apply to the moist dentine surface using a disposable brush.

Rub primer into the dentine for 30 s.

Remove excess primer with a brush and thoroughly dry the cavity with oil-free air.

Attention: Do not rinse with water!

5. Application of the Adhesive

Repeat this step a second time when using Rebuilda DC.

Dispense adhesive onto a mixing pallet and apply evenly to all prepared dentine and enamel surfaces using a new disposable brush. Thoroughly brush the adhesive into the primer for at least 15 s. Afterwards spread the adhesive gently with a faint air jet and polymerize with light for 20 s. To light-cure this material conventional polymerization devices are suited. The light output should be a minimum of 500 mW/cm² on halogen polymerization devices and 300 mW/cm² on LED devices. Do not remove or contaminate the unpolymerized inhibition layer as it is essential for bonding to the composite.

In order to produce a perfect bond between the hydrophilic dentine and hydrophobic composite, rubbing the primer into the smear layer (see item 4) and thoroughly brushing the adhesive into the primer/smear layer are essential.

6. Composite filling

Apply composite filling material and adapt it to the cavity walls. Build up composite filling layer by layer, polymerize and finish.

B) Direct composite restoration with total-etch-technique

Steps 1. and 2., see A) 1. and 2.

3. Etch the entire cavity with **Vococid** (etching gel); start with the enamel margins (enamel for ca. 30 s, dentine max. 15 s) remove the etching gel by suction, rinse for approx. 20 s and remove residual humidity (air syringe), but do not dry.

Continue 4. - 6., see A) 4. to 6.

C) Cementation of crowns, bridges, inlays, onlays and veneers

1. Preparation of the surfaces

Prepare tooth surfaces as usual. Mechanically remove remnants of cement, temporary luting materials as well as protective varnishes from the tooth surfaces to be treated, as otherwise **Solobond Plus** will not adhere. Prepare interior surfaces of the preparations to be bonded (sand-blasting, silanisè, etc.) according to the manufacturer's instructions for use. Afterwards avoid any contamination of the surfaces.

2. Etching of enamel margins, application of primer and adhesive

According to A) 1. - A) 5.

In case of tightly fitting preparations, make sure that the layers of the adhesive to be polymerized are not too thick in the area of the tooth surface.

3. Cementation

Insertion can be done either with chemically curing or light-curing composite luting cement. Apply cement according to the manufacturer's instructions for use. Insert and fix restorative.

Warnings/precautionary measures:

Phenolic substances (e.g. eugenol) affect polymerization and should not be used with **Solobond Plus**. Use a eugenol-free cement for luting temporaries.

Solobond Plus contains BIS-GMA, hydroethylmethacrylate, BHT, TEGDMA, acetone, sodium fluoride. Do not use in patients with known allergies to these ingredients. In case of linings containing composite materials or pulp capping materials, apply **Solobond Plus** intensively only to dentine areas (materials can be dissolved). Avoid contact of **Solobond Plus** with the oral mucosa. In case of contact of **Solobond Plus** with the eyes, rinse with plenty of water and consult a physician.

Information on storage and application:

Tightly close **Solobond Plus** container after use. Store at temperatures between 4°C - 23°C (storage at temperatures between 5°C - 8°C will extend the shelf life of the product). Do not expose **Solobond Plus** to direct sunlight or operating light. Do not use **Solobond Plus** after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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