VOCO

VOCO Ionofil Molar AC/Quick

(GB)

Instructions for use

in accordance with EN 29917/ISO 9917

(CE) 0482

VOCO Ionofil Molar AC/Quick is a radiopaque glass polyalkenoate (ionomer) cement for fillings in application capsules (AC) for direct application. **VOCO Ionofil Molar AC/Quick** adheres chemically to enamel and dentine and continuously releases fluoride ions.

VOCO Ionofil Molar AC is available in shades A1, A2, A3 and B3. **VOCO Ionofil Molar AC Quick** is available in shades A1, A2 and A3. From every application capsule 0,15 ml cement can be extruded.

Fields of application:


- restoration of V-shaped defects and cervical enamel or root erosions including class V cavities
- restoration of deciduous teeth (especially class I)
- fillings of class II cavities (not occlusion-bearing) by means of tunnel preparation
- build-up fillings and linings, core build-up
- retrograde root filling
- temporary long-term treatment of class I and II cavities

Application:**Preparation:**

The adhesion of **VOCO Ionofil Molar AC/Quick** to enamel and dentine prevents the formation of marginal gaps and creates retentions. Undercuts give additional stability to the filling. The layer thickness of fillings should not be less than 1 mm. Avoid feather-shaped margins. In case of doubt a conventional preparation should be preferred.

Clean the cavity (e. g. 3 % hydrogen peroxide), thoroughly rinse with water and dry (warm air). Do not overdry dentine. If the remaining dentine is thicker than 1,5 mm no additional lining is necessary. To avoid pulp irritations which may occur in case of insufficient dentine thickness, a covering of areas in proximity of the pulp with calcium hydroxide is recommended.

Before application, the capsule has to be activated, mixed and opened.

**Activation:**

For activation the liquid from the sachet underneath the clamp has to be completely pressed into the inner chamber of the capsule. This is done with the help of a capsule activator. Place the capsule in the activator so that the application cannula protrudes from it and the clamp faces the moveable part of the activator. **Hold the lever for 2 s in the compressed position.**

Mixing:

Fix the activated capsule in a capsule mixer (4.000 - 4.500 oscillations per min) and mix for approximately 10 s. Immediately remove the capsule from the mixer. The use of capsule mixers of other manufacturers may cause shorter or longer mixing times due to different mixing frequencies.

Opening:

Turn the twist cap with application cannula by a half turn to the right (see arrow on the twist cap) until you notice a stop.

Application:

Insert the opened capsule into an applicator. By turning the capsule in the applicator the cannula can be directed into the desired direction. Apply the material directly into the cavity by pressing the lever of the applicator within 30 s after mixing.

VOCO Ionofil Molar AC:

Mixing time	working time	setting time
10 s	2 min	4 min

VOCO Ionofil Molar AC Quick:

Mixing time	working time	setting time
10 s	1,5 min	2,5 min

After mixing, **VOCO Ionofil Molar AC/Quick** has a working time of approximately 1,5 min (room temperature). During application of a filling the working field must be kept dry. The use of a rubber dam is recommended.

Glass ionomer cements are susceptible to humidity (breath, saliva) or drying during the setting time, especially right after application. If possible, cover immediately with a strip or matrix or apply a protective varnish (see instructions for use **Final Varnish LC**).

Finishing and polishing (flexible polishing discs) of fillings should be conducted no earlier than 7 min (**VOCO Ionofil Molar AC**) respectively 2,5 min (**VOCO Ionofil Molar AC Quick**), after application of the filing. Recommendation: Re-apply protective varnish afterwards.

The best results of fine finishing will be achieved approximately 24 h after the filing has been inserted.

Additional indications/precautions:

VOCO Ionofil Molar AC/Quick should only be mixed immediately before application in order to have the full working time available for inserting and shaping and to avoid cracks from working the material too late.

Please note that the AC capsule has to be mixed right after its activation.

VOCO Ionofil Molar AC/Quick contains polyacrylic acid, fluoroallic acid and parabene, but does not use in cases of known sensitivities or allergies.

Storage

Store **VOCO Ionofil Molar AC/Quick** mixing capsules in a dry area at temperatures between 15°C - 23°C. Do not remove foil until ready to use. Do not use after the expiry date printed on the package.

Tightly close **Final Varnish LC** container immediately after use. Store at room temperature and avoid exposure to light.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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