Instructions for Use

V-Posil

Vinyl Polysiloxane Impression Material

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Instructions for Use

Product description:

V-Posil impression materials are addition-curing, elastomeric materials with high hydrophilic properties, high tear strength, high dimensional accuracy and a high resistance to permanent deformation.

The **V-Posil** family consists of different viscosities: putty, heavy-bodied, medium-bodied, light-bodied, x-light-bodied,

available in an assortment of delivery systems:

standard 1:1 (50 ml automix cartridges) and 5:1 (380 ml cartridges) for use in most commercially available automatic dispensing and mixing machines (e.g. **VOCO Dynamic Dispenser**, Sympress) and traditional 1:1 putty jars.

Indications/Techniques:

V-Posil Putty Fast is used as preliminary material for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- Two-step putty-wash impression technique using a foil (plastic putty spacer)
- One-step putty impression technique for forming functional peripheries

V-Posil Heavy Soft Fast is to be used as heavy bodied material for:

- One-step impression technique (simultaneous technique) using dual viscosities
- Two-step impression technique using dual viscosities
- Functional impressions

V-Posil Mono Fast is to be used as a medium bodied tray or syringeable impression material for:

- Taking impressions over fixed/removable restorations and implants (i.e., transferring impression posts and bridge components)
- Functional impressions
- Fabricating crown and bridgework or inlays
- Fabricating full or partial dentures
- Reline impressions
- Transferring root posts when fabricating posts and cores indirectly
- Multi tray technique
- Use in the simultaneous mixing technique as well as the putty-wash technique

V-Posil Light Fast and **V-Posil X-Light Fast** are to be used as syringeable impression materials for:

- Two-step putty-wash impression technique
- One-step putty-wash impression technique
- One-step impression technique using a foil (plastic putty spacer)
- One-step impression technique (simultaneous technique) using dual viscosities
- Reline impressions
- Fabricating full or partial dentures

Warnings:

Do not use **V-Posil** impression materials as a temporary reliner.

Do not use **V-Posil** impression materials with condensation-curing silicones, polyether, Vinylsiloxanether® or polysulfide materials.

V-Posil Putty Fast and **V-Posil Heavy Soft Fast** impression materials are not suitable for detailed impressions when used alone.

Precautionary measures:

Do not use after expiration date.

Do not leave any residual impression material in the sulcus or oral cavity.

Do not swallow impression material! If swallowed: consult a medical doctor if problems arise or persist.

Avoid contact with eyes. If accidental contact with the eyes occurs, rinse immediately and thoroughly with an eye wash or water.

Consult a medical doctor if problems arise or persist.

These products do not normally cause allergic reactions. However, for sensitive persons, an allergic reaction to the products cannot be ruled out.

Use of products containing active sulfur, aluminium chloride or nitrogen compounds (retraction cords containing ferric sulfate, polysulfide impression materials, etc.) in conjunction with this product will interfere with the setting reaction of the vinyl polysiloxane material. Use of the materials requires the area to be rinsed thoroughly to remove all residue, before an impression is taken. Do not use latex gloves.

Do not interchange the Base and Catalyst lids or scoops for hand-mixed jar putty.

In order to ensure an optimal adhesion of the two materials during the two-step putty-wash impression technique, both impression steps must take place one immediately after the other.

When taking impressions of areas with severe undercuts and wide interdental spaces, use standard block-out techniques.

When using a custom impression tray, ensure that sufficient space remains between the side of the tray and the teeth/jaw.

Do not allow the material to enter the sewer or water system to avoid environmental contamination.

Avoid contact with clothing, since the material cannot be removed by dry cleaning.

Note:

To ensure optimal impressions, the temperature of the material should not deviate from 23 °C (74 °F) before applying. Otherwise, working and setting times will be affected.

Dry prepared teeth before impression taking.

The warranted product characteristics can only be achieved when using the original Sulzer Mixpac mixing tips as approved by **VOCO**. All the mixing tips for use with our respective products are available from **VOCO**.

Impression tray: Preparation and adhesives

In principle, all common impression trays can be used if a relevant dynamic pressure is guaranteed. When retention is not sufficient, ensure a strong bond to the impression material, by brushing the impression tray with a thin film of **V-Posil Adhesive** for addition-curing silicone prior to loading the tray with impression material. Allow to dry per manufacturer's instructions.

Instructions for use:**Jar material for hand-mix putty preparation**

Only use materials with the same lot number.

Do not interchange the Base and Catalyst lids or scoops for hand-mixed jar putty.

Use the scoops to dispense equal amounts of Base and Catalyst.

The different color scoops must only be used for the materials with corresponding colors.

Close the containers carefully after use, and ensure that the lids are not interchanged.

Contamination of Base and Catalyst in the container renders the material unusable.

Knead the Base and Catalyst for 45 seconds until the material is homogeneously colored.

Should gloves be worn, test them for compatibility prior to mixing with a sample of the material to be mixed. Certain types of gloves (such as latex) can prevent the polymerization.

Using vinyl gloves is recommended.

Prior to taking the wash impression, the preliminary impression must be cleaned with water, dried and trimmed using standard methods.

Instructions for the safe handling of the VOCO Dynamic Dispenser system (380 ml cartridge system; 5:1) in dispensing device Sympress I or Sympress II

Prior to first use please unscrew the blue plug counter-clockwise and pull it out (figure 1, page 82). Check the outlet opening. Then screw on the dynamic mixer type 30 (figure 2, page 82). Ensure that the positioning aid is on the top.

Please use only the dispensing and mixing components compatible with the **VOCO** 380 ml cartridge system for dispensing and mixing impression material.

The **VOCO** 380 ml cartridge is heavy and could be damaged if dropped. Damaged cartridges must not be used under any circumstances.

When using the material in a Sympress II dispensing device, please activate the speed setting indicated on the product label. Failure to comply with this instruction may modify the curing characteristics of the impression materials. If you are using a unit other than the Sympress mixing device, please follow the instructions for the automatic dispensing and mixing system you are using. When changing to a different mixing device, always discard at least 3 cm of mixed material.

Using the VOCO Dynamic Dispenser system (380 ml cartridge system; 5:1) in dispensing device Sympress I or Sympress II

- By pressing the plunger return knob, allow the plunger to return to its original position automatically (while the cover is closed; figure 3, page 82). Next, open the cover (figure 4, page 82).
- Insert the **VOCO** 380 ml cartridge with the pre-assembled mixer in the Sympress I or Sympress II dispensing and mixing unit (figure 5, page 82). (If you are working with a previously used cartridge, the filled mixer will already be located on the cartridge as a cap).
- If you are using a Sympress II dispensing device, please activate the speed setting indicated on the product label. A light signal will confirm your selection.
- Close the device cover (figure 6, page 82).

Using the cartridge for the first time

- Dispense the impression material (ensure that the dynamic mixer type 30 is wound into place fully) until the base material and catalyst are mixed homogeneously (see coloring of the mixed impression material); discard at least the first 3 cm when dispensing (figure 7, page 82).
- At least 3 cm of the mixed material must be discarded each time the plunger is restarted or when using the cartridge in different mixing units.
- Load the impression tray with the required amount of material (figure 8, page 82) by pressing the appropriate feed key. Hold the impression tray at an angle when loading and press against it lightly. Leave the mixer in the material when dispensing to avoid air bubbles.
- Release the feed key as soon as the required amount of material has been dispensed. The filled dynamic mixer should remain on the cartridge as a seal.
- Remove the used mixer before the next application by turning it to the left (figure 9, page 82), visually check the outlet opening for blockage and attach a new dynamic mixer (figure 9, page 82).
- Remove the completely empty jumbo cartridge from the dispensing and mixing unit and dispose of it according to regional regulations.

Further use of a cartridge already in use

- Wind a new dynamic mixer onto the same mixing unit and begin the new application (material no longer has to be discarded). Proceed with the application as described above.

Instructions for use: Cartridge material for use in VOCO Dispenser Typ 2

Insert the cartridge into the **VOCO Dispenser Typ 2**. Ensure that the notches on the cartridge base are pointing down. The clasp will not close if the cartridge is not properly inserted into the dispenser.

Remove the cartridge cap. The cap can be replaced after initial use.

Advance plunger into cartridge and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a mixing tip on the cartridge and turn cap 90 degrees in a clockwise direction to lock in place.

Inject the required amount of material, e.g., directly into the tray or the preliminary impression. To apply the material around the prepared teeth use the cartridge with an intraoral tip type 3/type 6 on the mixer.

Replace the plug after use. Before using the cartridge again, remove the cartridge cap or the filled mixing tip and check the cartridge outlet openings to ensure that no polymerized materials have sealed the openings. Clear any clogged openings and dispense a small amount of impression material until equal amounts are extruded at the same rate.

Install a new mixing tip and continue as described above.

Storage:

Store in a dry place at room temperature; do not expose to sunlight and avoid exposure to extreme temperature fluctuations.

Disinfection:

Impressions can be disinfected using for example, a 2 % acidic glutaraldehyde solution. Use disinfection solutions specific for impression materials.

See manufacturer instructions.

Electroplating:

Impressions may be silver-plated with an alkaline-plating solution only.

Model preparation:

For model fabrication, the impression can be poured immediately after disinfection. Impressions can also be poured within two weeks with standard dental stone (class III - V).

Note:

Selected **VOCO** materials are available in certain markets only.

For professional use only. Caution: Federal (U.S.A.) Law restricts this device to sale by or on the order of a dentist.

GB

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

USA

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. VOCO recognizes its responsibility to replace products if proven to be defective. VOCO does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith. Descriptions and data constitute no warranty of attributes and are not binding.

CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

For questions or comments, please call 1-888-658-2584.

**Keep this material out of reach of children.
For dental use only.**

An explanation of the symbols used in labeling can be found at www.voco.com/us/symbols

Technical Data:

Product name	Mixing-ratio	Working time at 23 °C / 74 °F ≤	Working time at 35 °C / 95 °F ≤	Intraoral setting time at 35 °C / 95 °F ≥	Total setting time* ≥	Shore A hardness (approx.)	Linear dimensional change (approx.) %	Elastic Recovery test (approx.) %	Strain in compression (approx.) %
V-Posil Putty Fast	1:1	2 minutes		2 minutes	4 minutes	66	-0.25	98.5	3.0
V-Posil Heavy Soft Fast	5:1	2 minutes		2 minutes	4 minutes	55	-0.25	99.3	4.0
V-Posil Heavy Soft Fast	1:1	2 minutes		2 minutes	4 minutes	55	-0.25	99.3	4.0
V-Posil Mono Fast	5:1	2 minutes	1 minute	2 minutes	4 minutes	58	-0.25	99.5	4.0
V-Posil Mono Fast	1:1	2 minutes	1 minute	2 minutes	4 minutes	58	-0.25	99.5	4.0
V-Posil Light Fast	1:1	2 minutes	1 minute	2 minutes	4 minutes	47	-0.20	99.6	4.0
V-Posil X-Light Fast	1:1	2 minutes	1 minute	2 minutes	4 minutes	47	-0.20	99.6	4.5

*Total setting time (removal time from mouth) from start of mix

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putty, heavy-bodied, medium-bodied, light-bodied, x-light-bodied,

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VOCO GmbH