

# SCIENTIFIC INFORMATION

## Arabesk Flow – Technical Profile

Filler content	64.0 w/w %	DIN 51081
Polymerization shrinkage	3.47%	dilatometric
Compressive strength	330 MPa	analogous ISO 9917
3-Point flexural strength	113 MPa	ISO 4049
Modulus of elasticity	5800 MPa	ISO 4049
Radiopacity	200 %Al	ISO 4049
Colour stability	no discolouration	ISO 4049
Stability under daylight	meets standard	ISO 4049

The stated measurement values do not represent a specified value within the scope of the VOCO GmbH quality assurance system.