

## Ionoseal – Technical Profile

Depth of curing 20 s / 40 s / 60 s (> 500 mW / cm <sup>2</sup> )	5 mm / 10 mm / 10 mm	ISO 4049
Compressive strength	226 MPa	ISO 4049
Diametral tensile strength	42 MPa	
3-Point flexural strength	95 MPa	ISO 4049
Flexural modulus	3970 MPa	
Diffusion of fluoride	10 ppm/month	

The stated measurement values do not represent a specified value within the scope of the VOCO GmbH quality assurance system.